
PRINCIPIO DEL TEST

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
Gerhardt W., Waldenström, J. (1979). Clin.Chem., 25, 1274-1280.

PRINCIPLE

CK-MB is composed of two kind of subunits: M and B. A specific antibody inhibits the M subunit of the CK-MB and CK-MM not affecting the activity of the B subunit of the CK-MB and CK-BB isoenzymes. The remaining CK-MB activity, corresponding to the CK-B, is determined by the CK-NAC activated method.

The phosphate released from phosphocreatine by CK-B binds to ADP forming ATP, that reacts with glucose by the action of hexokinase producing glucose-6-phosphate and ADP.

The G-6-PDH degrades this product to phosphogluconate at the same time that NADP⁺ reduces to NADPH. This last reaction is measured by the variation of the Abs at 340 nm. In optimal conditions, the ΔAbs/min is proportional to the creatinase activity in the sample.

DIAGNOSTIC USE

The CK-MB increases take place in myopathic diseases and inflammatory diseases of the heart. In myocardial infarction, CK-MB activity increases rapidly. When serial determinations are performed it is observed that the maximum point is quickly reached with an also very quick decrease.

Single test result can not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 20 x 2.5 mL (Ref. 99 91 11). Contents

A. 20 x 2.5 mL Substrate / antibody	Ref. 99 60 52
B. 1 x 50 mL Buffer solution	Ref. 99 43 35
C. 1 x 2 mL Control	Ref. 99 30 56

Concentration is stated on the label

PREPARATION OF WORKING REAGENT AND CONTROLS

Reagent: Rehydrate 1 vial of freeze-dried substrate/antibody with the volume of buffer solution stated on the label. Mix gently until completely dissolved.

Control: Reconstitute the contents of the vial with the volume of deionised water stated on the label. Leave to stand for 15 min and dissolve completely by mixing gently.

WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:

Imidazole/Acetic acid buffer pH 6.8	100 mM
Glucose	20 mM
Creatine phosphate	30 mM
Mg ²⁺	10 mM
ADP	2 mM
AMP	5 mM
di-adenosine pentaphosphate	10 μM
N-acetyl cysteine	20 mM
NADP ⁺	2 mM
Glucose-6-phosphate dehydrogenase	≥ 1,500 U/L
Hexokinase	≥ 2,500 U/L
Antibody to CK-M	≥ 2,000 U/L
Stabilizers and preservatives	

Control: Pool of human sera with a known CK-MB enzymatic activity, added with stabilizers and preservatives

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The rehydrated reagent is stable for 12 days at 2-8°C and for 3 days at room temperature (≤ 25°C). The reconstituted control is stable for 8 hours at room temperature and 5 days at a 2-8°C. If deep frozen, at -20°C it will remain stable for 4 weeks. Freeze and thaw only once.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank ≥ 0.800.

ADDITIONAL EQUIPMENT

General laboratory equipment.
Spectrophotometer; automated analyzer or photometer with thermostated cuvette, 1 cm light path.

SAMPLE

Serum or plasma, with heparin.
When the sample is stored at 2°-8° C, the CK - MB activity remains stable for one week..

CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.
The safety statements are on the label. It is advisable to look at the SDS before using the reagent.
The disposal of the residues has to be made according to local regulations.

QUALITY CONTROL

It is advisable to include the CK-MB control (Ref. 99 30 56) along with other sera, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) in each test series for results verification. Each particular laboratory should establish its own control program.

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

Prior to assay the CK-MB, the CK total activity has to be determined (QCA Ref. 99 44 10) or CK-NAC Liquid (Ref. 99 05 24 or 99 79 74). If the CK total activity ≥ 1000 U/L, dilute the sample 1/10 with saline (NaCl 0.9%), and then proceed with the CK-MB test.

Technique

Bring reagents and the analyzer to the desired assay temperature (25°, 30°, 37°C).

Technique	mL
Working reagent	1.0
Incubate 2-3 min at the desired assay temperature	
Sample	0.04

Mix well and let stand 10 min. at the assay temperature. Then add the mixture to a cuvette and read Abs₁, and then read Abs₂ exactly 5 min. later. Determine: ΔAbs = Abs₂ - Abs₁

Reading

Wavelength: 365 nm; 340 nm; 334 nm
Blank: water
Cuvette: thermostated 1 cm light path

CALCULATIONS

Using the stated formula for factor and U/L calculation:

$$\Delta \text{Abs} \times \frac{\text{Vt} \times 10^6}{\text{C} \times \text{l} \times \text{Vs} \times \text{t}} = \text{U/L}$$

Vt: Total volume of reaction mixture;
Vs: Sample volume
l: Cuvette light path
C: Extinction coefficient of NADPH:
365 nm: 3.53 x 10³
340 nm: 6.31 x 10³
334 nm: 6.17 x 10³

t: Time of reading

To calculate the U/L with the values of ΔAbs

	CK-B	CK-MB
365 nm	1486 x ΔAbs	2972 x ΔAbs
340 nm	825 x ΔAbs	1651 x ΔAbs
334 nm	841 x ΔAbs	1683 x ΔAbs

As the activity tested is the one due to CK-B subunit, to retrieve the value of the total activity of the isoenzyme CK-MB the U/L obtained are to be multiplied by 2. This is why the factor used to calculate CK-MB activity is double.

REFERENCE VALUES

There will be suspicion of myocardial infarction, provided the following three conditions are given:

	25°C	30°C	37°C
CK men	> 80 U/L	> 130 U/L	> 190 U/L
CK women	> 70 U/L	> 110 U/L	> 170 U/L
CK - MB	> 10 U/L	> 16 U/L	> 25 U/L

$$\text{CK-MB index} = \frac{\text{CK-MB activity}}{\text{CK total activity}} \times 100: 6 - 25\%$$

The stated values are for guidance. Each particular laboratory should establish its own normal range, using its own instrumentation, blood collection methods and test procedures

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using an automated method at 37°C and 340nm.

Sensitivity, as detection limit: 4 U/L

Linearity: For values higher than 600 U/L, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.

Repetitivity, as Coefficient of Variation: 2.19%

Reproducibility, as Coefficient of Variation: 2.52%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

INTERFERENCES

Highly haemolysed sera will interfere with the reaction. Glucose up to 600 mg/dL, ascorbic acid up to 40 mg/dL and haemoglobin up to 500 mg/dL, do not interfere the assay.

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GB - CK-NAC / CK-MB, BUFFER

Danger

Hazard: May damage the unborn child.
Precautionary: Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves/protective clothing/eye protection/face protection. IF exposed or concerned: Get medical advice/attention. Store locked up. Dispose the contents/container according to point 13 of the Safety Data Sheet.

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