# LDH

## MÉTODO SFBC

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# LDH

## SFBC METHOD

For "in vitro" determination of lactate dehydrogenase in serum or plasma

### PRINCIPLE

The lactate dehydrogenase catalyzes the reduction of pyruvate with the aid of NADH as cofactor, which is oxidized to NAD<sup>+</sup>, giving an absorbance change. In optimum reaction conditions the  $\Delta\text{Abs}/\text{min}$  is directly related to LDH concentration in the sample.

### DIAGNOSTIC USE

Increased levels of LDH in serum are found in patients suffering from hepatitis, infectious mononucleosis, malignant tumours, pancreatitis, cirrhosis, or myocardial infarction, due to release of enzyme from the tissue. Values also increase in cases of muscular dystrophy or anemias. Values below the normal range are not clinically significant.

Single test result can not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

Kit 20 x 3 mL. (Ref. 99 82 87). Contents:

- A. 20 vials of freeze-dried NADH
- B. 1 x 60 mL Buffered substrate

Ref. 99 17 67  
Ref. 99 91 25

### WORKING REAGENT PREPARATION

Rehydrate one vial of NADH with the volume of buffered substrate stated on the label and mix gently.

### WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:

Tris-HCl buffer pH 7.2	80 mM
Sodium pyruvate	1.6 mM
NaCl	200 mM
NADH	0.20 mM
Stabilizers and preservatives	

### STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The reconstituted reagent is stable for 21 days at 2-8°C and for 5 days at room temperature ( $\leq 25^\circ\text{C}$ ), when protected from the sunlight.

### Signs of reagent deterioration:

Reagent turbid or with visible particles. Working reagent blank:  $\leq 1.0$

### ADDITIONAL EQUIPMENT

General laboratory equipment.  
Spectrophotometer; automated analyzer or photometer with thermostated cuvette, 1cm light path.

### CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.  
The safety statements are on the label. It is advisable to look at the SDS before using the reagent.  
The disposal of the residues has to be made according to local regulations.

### SAMPLE

Serum, heparinized plasma. Samples free from hemolysis should be used. The enzyme in serum is stable for 2 days at 2-8°C. Frozen samples are rapidly inactivated.

### QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

### AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

### INTERFERENCES

It is recommended to use heparin as anticoagulant. Other anticoagulants, like citrate, oxalate, or fluoride may interfere with the test. Hemolyzed samples will give false results.  
When assaying high activity samples, a very low initial absorbance can be found, which is mainly due to the rapid consumption of the NADH at the early stage of the reaction. In such a case, dilute the sample with saline (NaCl 0.9 %) and assay it once again.

### PROCEDURE

The method described is the one proposed by the French Society of Clinical Biology. Bring the working reagent and the instrument to the working temperature (25, 30 or 37°C).

Technique (25°C/30°C)	Macro (mL)	Semimicro ( $\mu\text{L}$ )
Working reagent	3.0	500
Incubate 2-3 min. at 25 or 30°C		
Sample	0.1	20

Technique (37°C)	Macro (mL)	Semimicro ( $\mu\text{L}$ )
Working reagent	3.0	500
Incubate 2-3 min. at 37°C		
Sample	0.05	10

Mix and start the stopwatch. Transfer to the measuring cuvette.  
Read the absorbance after 1, 2, 3 and 4 min. Determine the average value of  $\Delta\text{Abs} / \text{min}$ .

### Reading

Wavelength: 334 nm; 365 nm; 340 nm  
Blank: Water  
Cuvette: Thermostated 1 cm light-path

### CALCULATIONS

The formula indicated is used to obtain the factor to calculate the U/L

$$\Delta\text{Abs}/\text{min} \times \frac{Vt \times 10^3}{\epsilon \times l \times V_s} = \text{U/L}$$

Where:

Vt: Total volume of reaction mixture;  
Vs: Sample volume  
l: Cuvette light path  
 $\epsilon$ : Extinction coefficient of NADH:

365 nm: $3.40 \times 10^3$
340 nm: $6.31 \times 10^3$
334 nm: $6.17 \times 10^3$

Determine the average value of  $\Delta\text{Abs} / \text{min}$ .

$$\text{U/L} = \Delta\text{Abs}/\text{min} \times \text{Factor}$$

Factors

		25°C - 30°C	37°C
Macro	334 nm	5016	9871
	340 nm	4921	9682
	365 nm	9118	17941
Semimicro	334 nm	4207	8252
	340 nm	4127	8095
	365 nm	7647	15000

### REFERENCE VALUES

Temperature	Adults
25 °C	120 - 240 U/L
30 °C	160 - 320 U/L
37 °C	225 - 450 U/L

The stated values are for guidance. Each particular laboratory should establish its own normal range, using its own instrumentation, blood collection methods and test procedures

### PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using a manual method at 37°C and 340nm.

Sensitivity, as detection limit: 10 U/mL.

Linearity: Up to 2300 U/L. For higher values, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.

Accuracy: 98.6%

Repetitivity, as Coefficient of Variation: 2.0%

Reproducibility, as Coefficient of Variation: 3.4%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

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