# COLESTEROL - HDL

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Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

Zapico E. y Ordóñez J. (2002), Clin Invest Ar ter ioscl;14(5):272-6.

# HDL - CHOLESTEROL

## DEXTRAN SULPHATE - Mg (II) METHOD

For "in vitro" determination of HDL - cholesterol in serum

### PRINCIPLE

LDL and VLDL (low and very low density lipoproteins) are precipitated from serum by the action of a sulfated polysaccharide, in the presence of divalent cations. Then, the cholesterol bound to high density lipoproteins (HDL- cholesterol) present in the supernatant, is determined.

### DIAGNOSTIC USE

The fraction of the cholesterol bound to high density lipoprotein is an indicator of the risk of coronary heart disease. High levels of HDL-cholesterol appear to act as a protective factor, while low values are one of the major risk factors.

Determination of HDL- cholesterol, along with comprehensive study of the lipid profile of the patient, allows assessing the risk of coronary heart disease.

Low levels of HDL- cholesterol are found in cases of unbalanced diet, sedentary lifestyle, alcoholism or smoking.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

1 x 4 mL Precipitant solution  
Dropper for a minimum of 100 tests

Ref. 99 38 85

### WORKING REAGENT PREPARATION

Reagent is ready to use

### WORKING REAGENT COMPOSITION

Concentrations in the reagent solution are:

Dextran sulphate	10 g/L
Magnesium acetate	1 M
Stabilizers	

### STORAGE AND STABILITY

When kept at 2 - 8°C, the reagent will remain stable until the expiration date stated on the label.

### Signs of reagent deterioration:

Presence of particles or turbidity in the reagent.

### ADDITIONAL EQUIPMENT

General laboratory equipment.  
Spectrophotometer or photometer thermostable at 37°C.  
Centrifuge.

Reagent for cholesterol determination and standard  
It is advisable to use QCA reagents:

Kit 1 x 100 mL 99 52 82  
Kit 3 x 100 mL 99 52 80  
Kit 2 x 250 mL 99 50 12

Standard is included in the kits

### CAUTION

The safety statements are on the label.  
We advise to read MSDS before reagent handling. Waste products must be handled as per local regulations.

### SAMPLE

Serum. The separation of the HDL fraction has to be carried out as soon as possible; otherwise, it is recommended to freeze the sample at - 15°C. Stored in this way, it will remain stable for up to 1 week.

### QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series for determining cholesterol.

For the precipitation reaction control, it is recommended to use liquid or lyophilized serum in which the native state of lipoproteins is retained. We recommend using a "pool" of normal sera of known value in HDL-cholesterol.

Each particular laboratory should establish its own control program.

### AUTOANALYZERS

The precipitation reaction cannot be automated  
Adaptations for cholesterol determination to different autoanalyzers are available on request.

### PROCEDURE

#### 1. Precipitation reaction:

Sample 0.3 mL  
Precipitant solution 1 drop  
Mix and let stand for 15 min at room temperature (20 - 25°C).  
Centrifuge at 2000 x g (1500 -2300 rpm) / 15 min or  
Centrifuge 10000 x g (8000 - 12000 rpm) / 2 min  
Determine, in the supernatant, the concentration of cholesterol.

#### 2.- Determination of cholesterol (1)

Technique	BL (mL)	SA (mL)	ST (mL)
Supernatant	--	0.01	--
Standard	--	--	0.01
Working reagent	1.00	1.00	1.00

Mix well and let stand for 5 min at 37°C or 10 min at room temperature.

#### Reading

Wavelength: 546 nm; 505 nm  
Blank: the contents of BL  
Colour stability: 1 hour (when protected from direct light)

#### CALCULATIONS

$\frac{\text{Abs. SA}}{\text{Abs. ST}} \times 200 \times 1.13 = \text{mg HDL- cholesterol / dL}$

Where:

Abs SA: Sample Absortion  
Abs ST: Standard Absortion

S.I. Units: mg/dL x 0.0259 = mmol / L

**Note: When the precipitant reagent is added to the sample, the latter is diluted by a factor of 1.13. Therefore, the cholesterol value should be multiplied by such a factor. (1)With QCA cholesterol reagent**

#### REFERENCE VALUES

##### Coronary heart disease risk:

< 40 mg/dL High  
≥ 60 mg/dL Low

The concentrations vary considerably with age and sex.

#### Performance characteristics

The performance characteristics depend on the precipitant reagent, the reagent to determine cholesterol and the method used. These results have been obtained by manual method. The cholesterol in the supernatant has been determinate with the reagent cholesterol Liquid from QCA.

Sensitivity, as detection limit: 3 mg/dL

Accuracy: 97.4%

Repetitivity, as CV%: 1.87%

Reproducibility, as CV%: 2.22%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

#### INTERFERENCES

Neither aged nor haemolysed samples shall be assayed. Bilirubin, at concentrations higher than 9 mg/dL, interferes in the precipitation reaction.

Sera with triglycerides concentration higher than 350 mg/dL will be diluted 1/2 in saline (NaCl 0.9%) prior to precipitation.

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Albers,J.J., Warmick,G.R., Cheng,M.C. (1978). Lipids, 13, 926 - 932.

Benzie,I. (1979). Med. Lab. Sci., 36, 289 - 291.

Wieland, H. Seidel, D. (1981). Ärtztl. Lab.,27, 141 - 154.

A policy statement of the European Atherosclerosis Society, European Heart Journal 8,(1987) 77 - 88.

Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

National Cholesterol Educacion Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

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