GLUCOSA LÍQUIDA

MÉTODO GOD – POD

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Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACCPress (2000).

Trinder, p. (1969). Ann Clin. Chem. 6, 24 - 27.

GLUCOSE LIQUID


GOD – POD METHOD

For "in vitro" determination of Glucose in serum, plasma or CSF

PRINCIPLE

The oxidation of glucose to gluconic acid is catalyzed by glucose oxidase producing hydrogen peroxide.

The hydrogen peroxide reacts with 4-aminoantipyrine and p-hydroxybenzoic acid in the presence of peroxidase to give a quinone derivative, whose coloration is proportional to the glucose concentration in the sample.

DIAGNOSTIC USE

The determination of glucose in serum or urine is used for the evaluation of disorders of the metabolism of carbohydrates.

Glucose is the major source of energy for body cells. Insulin produced in the pancreatic cells, facilitates the entry of glucose into cells of tissues.

Increased blood glucose is associated with a decrease in insulin activity or with its deficiency.

In serum or plasma it is found elevated glucose values mainly in patients with diabetes mellitus but also acute pancreatitis, Cushing syndrome, acromegaly and gigantism.

Hypoglycemia can occur in response to fasting, or may be due to drugs, poisons or inborn errors of metabolism.

The presence of glucose in urine with the individual not suffering from diabetes is usually a sign of disease in the renal tubules.

Glucose determination in CSF is primarily of interest in case of bacterial meningitis, when its concentration is low or not detected.

Single test result cannot be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 1 x 100 mL. (Ref. 99 82 25) Contents:

A. 1 x 100 mL Reagent Ref. 99 82 84
B. 1 x 5 mL Standard Ref. 99 02 93

Kit 3 x 100 mL. (Ref. 99 82 82) Contents:

A. 3 x 100 mL Reagent Ref. 99 82 84
B. 1 x 5 mL Standard Ref. 99 02 93

Kit 4 x 250 mL. (Ref. 99 86 60) Contents:

A. 4 x 250 mL Reagent Ref. 99 01 68
B. 1 x 5 mL Standard Ref. 99 02 93

WORKING REAGENT PREPARATION

Reagent and standard are ready to use.

REAGENT COMPOSITION

Concentration in the reagent solution is:

Phosphate buffer pH 6.8	100 mM
p-hydroxybenzoic acid	39.5 mM
4-aminoantipyrine	0.8 mM
Phenol	4.5 mM
Glucose oxidase	≥ 18 kU/l
Peroxidase	≥ 1.1 kU/l
Preservatives and stabilizers	

Standard: Aqueous solution equivalent to 100 mg/dL (5.55 mmol/L).

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank >0.400.

ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer thermostatable at 37°C with a 505 nm filter. Cuvette: 1 cm light-path.

SAMPLE

Serum, plasma or CSF. Serum or plasma glucose (but not whole blood glucose, due to glycolytic processes) can be stored during 2 - 3 days, when refrigerated at 2-8°C.

The CSF must be clear and without debris. In these conditions the glucose is stable 48 hours at 2-8°C.

CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the SDS before using the reagent.

The disposal of the residues has to be made according to local regulations.

PROCEDURE

Bring the reagent and the analyzer to the working temperature.

Technique	BL mL	SA mL	ST mL
Standard	--	--	0.01
Sample	--	0.01	--
Working reagent	1.0	1.0	1.0

Mix well and incubate 5 - 10 min. at 37°C or 20-25 min. at 20 - 25°C.

Reading

Wavelength: 505 nm.

Blank: the contents of BL.

Colour stability: a minimum of 1 hour, when protected from direct sunlight.

CALCULATIONS

$\frac{\text{SA Abs.}}{\text{ST Abs.}} \times 100 = \text{mg glucose / dL}$

ST Abs.

Where:

SA Abs: Sample Absorption

ST Abs: Standard Absorption

S.I. Units

(mg/dL) x 0.0555 = mmol/L

REFERENCES VALUES

Serum, plasma (fasting):

Adult: 74 - 115 mg/dL (4.1-6.4 mmol/L)

Child: 60 - 100 mg/dL (3.3-5.6 mmol/L)

Neonate: 30 - 80 mg/dL (1.7-4.5 mmol/L)

Neonate premature: 20 - 60 mg/dL (1.1-3.3 mmol/L)

CSF:

Adult: 40 - 70 mg/dL (2.2-3.9 mmol/L)

Child: 60 - 80 mg/dL (3.3-4.5 mmol/L)

Urine: 1 - 15 mg/dL (0.1-0.8 mmol/L)

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using a manual method.

Sensitivity, as detection limit: 2.0 mg/dL

Linearity: Up to 500 mg/dL. For higher values, it is recommended to dilute the sample 1/2 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 2.

Accuracy: 98.9 %.

Repeatability, as Variation Coefficient: 0.79%

Reproducibility, as Variation Coefficient: 1.33%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request

INTERFERENCES

Haemoglobin, higher than 200 mg/dL; Bilirubin, higher than 20 mg/dL; Uric acid, higher than 20 mg/dL; Creatinine, higher than 15 mg/dL.

Interferences caused by the anticoagulants of current use such as Heparin, EDTA or oxalate have not been described.

QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

REFERENCES

Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, PA

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

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