CREATININA

MÉTODO DE JAFFE MODIFICADO

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Spierto, F.W., McNeil, M.L., Burtis, C.A. (1979), Clin. Biochem. 12, 18 - 21.

Clinical Diagnosis and Management, Todd-Sanford-Davidsohn, Ed. J.B. Henry 6ª Edition 1979. W.B. Saunders, Philadelphia, 262 - 263.

CREATININE

MODIFIED JAFFE METHOD

For "in vitro" determination of creatinine in serum, plasma or urine

PRINCIPLE

In alkaline medium, creatinine with picric acid form a colored compound, alkaline creatinine picrate which is determined photometrically. The color produced in the reaction is proportional to the concentration of creatinine in the sample under optimal assay conditions.

DIAGNOSTIC USE

Creatinine in serum increases in case of renal failure, acute or chronic, hyperthyroidism, acromegaly and gigantism.

Increasing creatinine in urine occurs in diabetes mellitus, and infection. Also, exercise promotes increased excretion in the urine.

During pregnancy or decreased muscle mass, serum creatinine concentration is reduced, while creatinine in urine is decreased in renal failure, myopathy, anemia or hypothyroidism.

One single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

Reagents

Kit 2 x 100 mL. (Ref. 99 88 91). Contents:

- A. 1 x 100 mL Picric acid solution Ref. 99 90 06
- B. 1 x 100 mL Alkaline solution Ref. 99 56 35
- C. 1 x 5 mL Standard Ref. 99 83 99

Kit 4 x 250 mL. (Ref. 99 00 77). Contents:

- A. 2 x 250 mL Picric acid solution Ref. 99 00 82
- B. 2 x 250 mL Alkaline solution Ref. 99 01 08
- C. 1 x 5 mL Standard Ref. 99 83 99

WORKING REAGENT PREPARATION

Mix equal parts of each reagent (A and B) prior to assay.

WORKING REAGENT COMPOSITION

The reagent composition is as follows:

Picric acid	55 mM
Sodium carbonate	50 mM
NaOH	0.40 M
Preservatives and stabilizers	

Standard. Aqueous solution equivalent to 2 mg/dL, (176.8 µmol/L). Ready-to-use.

STORAGE AND STABILITY

The components of the kit, when stored at room temperature ($\leq 25^{\circ}\text{C}$), will remain stable until the expiration date stated on the label.

Working reagent is stable for up to 15 days at room temperature ($\leq 25^{\circ}\text{C}$). Store in the dark.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent.

Working reagent blank $\geq 0,300$.

ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer thermostatable at 37°C . Cuvette: 1 cm light-path

SAMPLE

Serum or plasma with heparin and urine. Serum or plasma samples are stable for 24 hours at $2 - 8^{\circ}\text{C}$.

To assay urinary creatinine, dilute the sample 1/20 in deionized water. Multiply the final result by 20.

CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the MSDS before using the reagent.

The disposal of the residues has to be made according to local regulations.

QUALITY CONTROL

It is recommended to include control sera, Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Abnormal (Abnormal Control Serum) (Ref. 99 46 85), in each test series for assuring the obtained values .

Each particular laboratory should establish its own quality control program and measures for avoiding results deviations

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

Bring the reagent to working temperature

Technique	BL	SA	ST
	mL	mL	mL
Standard	---	---	0.10
Sample	---	0.10	---
Working reagent	1.00	1.00	1.00

Mix and start the stopwatch. Transfer to the measuring cuvette.

Read the absorbances at 20 and 80 seconds

Reading

Wavelength: 546 nm; 510 nm

Blank: the contents of BL

CALCULATIONS

a) Creatinine concentration:

Determine the $\Delta\text{Abs.}$ for sample and standard.

$$\Delta\text{Abs.} = \text{Abs } 80 \text{ sec} - \text{Abs } 20 \text{ sec}$$

$$\frac{\Delta\text{Abs. SA}}{\Delta\text{Abs. ST}} \times 2 = \text{mg creatinine/dL}$$

b) To calculate the Clearance:

$$\frac{(\text{mg creatinine/dL urine}) \times \text{mL urine}/24 \text{ h.}}{(\text{mg creatinine/dL serum}) \times 1440} = \text{mL / min}$$

Where:

mg/dL urine and serum: values obtained in a)

mL urine/24h: excreted urine volume in 24h

1440: conversion factor (24 hours to minutes)

S.I. Units

$(\text{mg/dL}) \times 88.4 = \mu\text{mol/L}$

REFERENCES VALUES

	Serum, plasma	Urine (*)
Men:	0.6-1.1 mg/dL	21-26 mg/Kg/24 h.
Women:	0.5-0.9 mg/dL	6-22 mg/Kg/24 h.

Clearance

Men: 97-137 mL/min

Women: 88-128 mL/min

(*)The values are indicated as amount of creatinine excreted per Kg of weight in 24h.

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

PERFORMANCE CHARACTERISTICS

Performance of the reagent depends on the reagent itself and also depends on the method and the analyzer used.

The results indicated have been obtained using a manual method.

Sensitivity, as detection limit: 0.03 mg/dL

Linearity: Up to 15 mg/dL of Creatinine. For higher values, dilute the sample 1:1 in deionized water and assay once again. Multiply the final result by 2.

Accuracy: 97.4%

Repetitivity as Variation Coefficient: 1.72%

Reproducibility, as Variation Coefficient: 2.11%

Trueness: Results obtained with this reagent did not show sensible differences when compared with the reference reagent.

Details of the performance studies are available on request

INTERFERENCES

Lipemic and hemolyzed sera will interfere with the assay.

Several medicinal substances, eventually present in the sample, such as ascorbic acid or levodopa, can interfere with the final result.

At the described assay conditions, protein and glucose interferences are minimized.

REFERENCES

Blass, K.G., Thibert, R.J., Lam, L.K.(1974) Z. Clin. Chem. Clin. Biochem. 12, 336 - 343.

Spierto, F.W., McNeil, M.L., Burtis, C.A. (1979), Clin. Biochem. 12, 18 - 21.

Clinical Diagnosis and Management, Todd - Sanford - Davidsohn, Ed. J. B. Henry 6th Edition 1979. W. B. Saunders, Philadelphia, 262 - 263.

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GB - CREATININE (C/S DEPROTEINIZATION), PICRIC

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GB - CREATININE (JAFFE), ALKALINE SOLUTION

Warning

Hazard: Causes skin irritation. Causes serious eye irritation.

Precautionary: Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

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