# BILIRRUBINA

## MÉTODO JENDRASSIK – GROF

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# BILIRUBIN

## JENDRASSIK - GROF METHOD

For "in vitro" determination of direct and total bilirubin in serum or plasma

### PRINCIPLE

Total bilirubin is determined by reaction with diazotized sulfanilic acid, in the presence of caffeine, with the final production of an azopigment. The same reaction, but in the absence of caffeine, is used to measure direct bilirubin.

### DIAGNOSTIC USE

Bilirubin is a metabolic hemoglobin product which is transported to the liver, where is conjugated with glucuronic acid and it is excreted through bile. This "conjugated" bilirubin is called direct bilirubin, while unconjugated bilirubin is called indirect bilirubin. The total serum bilirubin is direct bilirubin plus indirect bilirubin.

Direct bilirubin increases in biliary obstruction, cirrhosis, and hepatitis. High level of indirect or total bilirubin may indicate hemolytic diseases, physiological neonatal jaundice, Gilbert's disease or fructose intolerance.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

Kit 1 x 245 mL. (Ref. 99 20 93). Contents:

A. 1 x 40 mL Sulphanilic acid solution.

B. 1 x 100 mL Caffeine solution.

C. 1 x 100 mL Tartrate solution.

D. 1 x 5 mL Sodium nitrite solution.

Ref. 99 29 04

Ref. 99 23 82

Ref. 99 29 46

Ref. 99 25 25

### WORKING REAGENT PREPARATION

All reagents are ready to use.

### REAGENTS COMPOSITION

The concentrations in the reagents solution are:

A. Sulphanilic acid 30 mM  
HCl 0.17 M

B. Caffeine 0.25 M  
Sodium benzoate 0.50 M

C. Potassium tartrate 0.80 M  
NaOH 1.85 M

D. Sodium nitrite 30 mM

### STORAGE AND STABILITY

When kept at room temperature ( $\leq 25^{\circ}\text{C}$ ), the components of the kit will remain stable until the expiration date stated on the label.

### Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank  $> 0,100$ .

### ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer; automated analyzer or photometer. Cuvette: 1 cm light path

### CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.

The safety statements are on the label. It is advisable to look at the MSDS before using the reagent.

The disposal of the residues has to be made according to legal local regulations.

### SAMPLE

Serum or plasma free from hemolysis. Serum bilirubin will decrease 50% in one hour if stored at room temperature ( $\leq 25^{\circ}\text{C}$ ) and at direct sunlight. Samples stored in the dark at  $2 - 8^{\circ}\text{C}$  will remain stable for up to 3 days.

### QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

### REFERENCE VALUES

Total bilirubin.	Adults: up to 1.0 mg/dL		
	Neonate:	Premature	Full term
	< 1 day	1.0 - 8.0 mg/dL	2.0 - 6.0 mg/dL
	< 2 days	6.0 - 12.0 mg/dL	6.0 - 10.0 mg/dL
	3-5 days	10.0 - 14.0 mg/dL	4.0 - 8.0 mg/dL
Direct bilirubin	Adults: up to 0.25 mg/dL		

The stated values are for guidance. Each particular laboratory should establish its own normal range.

### PROCEDURE

#### TOTAL BILIRUBIN

Technique	BL mL	SA mL
Sulphanilic acid	0.2	0.2
Sodium nitrite	--	1 gota
Caffeine	1.0	1.0
Sample	0.2	0.2
Mix and let stand 10 min. at room temperature ( $20 - 25^{\circ}\text{C}$ ).		
Tartrate	1.0	1.0
Mix and let stand for 5 min at room temperature ( $20 - 25^{\circ}\text{C}$ ).		

#### Reading

Wavelength: 578 nm.

Blank: BL contents

Cuvette: 1 cm light-path.

Color stability: a minimum of 1 hour.

#### DIRECT BILIRUBIN

Technique	BL	SA
Sulphanilic acid	0.2 mL	0.2 mL
Sodium nitrite	--	1 drop
Saline	2.0 mL	2.0 mL
Sample	0.2 mL	0.2 mL
Mix and let stand for 5 min at room temperature ( $20 - 25^{\circ}\text{C}$ ).		
<b>Read after exactly 5 min.</b>		

#### Reading

Wavelength: 546 nm.

Blank: BL contents

Cuvette: 1 cm light-path.

#### CALCULATIONS

Abs. SA x 10.8 = mg total bilirubin / dL.

Abs. SA x 14.4 = mg direct bilirubin / dL.

Where: Abs SA= Sample absorption

#### SI Units

(mg/dL) x 17.1 =  $\mu\text{mol/L}$ .

**NOTE:** In case of small-volume samples, please use 50  $\mu\text{L}$ , in place of the 0,2 mL as herein indicated, but keeping the same reagent volumes, and multiply the results by the following factors:

**F = 43.2** (total bilirubin)

**F = 57.6** (direct bilirubin)

#### PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. The indicated results have been obtained using a manual method.

Sensitivity, as detection limit: 0.1 mg/dL

Linearity: Up to 25 mg/dL. Concentrations higher than 15 mg/dL may produce an absorbance value (Abs  $> 1.500$ ) that can hardly be read in non-digital photometers. It is recommended in such a case, to carry out a 1/10 dilution of the sample with saline, (NaCl 0,9%) and then to multiply the final result by 10.

Accuracy: 96.3 %

Repetitivity, as Variation Coefficient: 1.64%

Reproducibility, as Variation Coefficient: 1.98%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

#### INTERFERENCES

Hemolyzed samples could give false results.

#### REFERENCES

Jendrassik, L., Grof, P.(1938). Biochem. Z., 297, 81-89.

Doumas, B. T., Perry, B. W., Sasse, E. A., Straumfjord Jr., J.V.(1973). Clin.Chem., 19, 984-993.

Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

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Tietz, N.W., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

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GB - BILIRUBIN (J-G METHOD), TARTRATE

#### Danger

Hazard: Causes severe skin burns and eye damage.  
Precautionary: Do not breathe dust/fume/gas/mist/vapours/spray. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Specific treatment.

Contains: caustic soda,sodium hydroxide

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