LDH LÍQUIDA

MÉTODO SFBC

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LDH LIQUID

SFBC METHOD

For "in vitro" determination of lactate dehydrogenase in serum or plasma

PRINCIPLE

The lactate dehydrogenase catalyzes the reduction of pyruvate with the aid of NADH as cofactor, which is oxidized to NAD⁺, giving an absorbance change. In optimum reaction conditions the $\Delta\text{Abs}/\text{min}$ is directly related to LDH concentration in the sample.

DIAGNOSTIC USE

Increased levels of LDH in serum are found in patients suffering from hepatitis, infectious mononucleosis, malignant tumours, pancreatitis, cirrhosis, or myocardial infarction, due to release of enzyme from the tissue. Values also increase in cases of muscular dystrophy or anemias. Values below the usual are not clinically significant.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

Kit 1 x 50 mL. Ref. 99 18 18. Contents:

A. 1 x 40 mL Buffer solution Ref. 99 13 74
B. 1 x 10 mL Buffered NADH Ref. 99 32 19

Kit 1 x 125 mL. Ref. 99 00 35. Contents:

A. 1 x 100 mL Buffer solution Ref. 99 00 40
B. 1 x 25 mL Buffered NADH Ref. 99 00 45

WORKING REAGENT PREPARATION

Reagents A and B are ready-to-use. If a monoreagent procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of A (buffer solution) + 1 part of B (buffered NADH).

WORKING REAGENT COMPOSITION

The concentrations in the reagent solution are:

Tris-HCl buffer pH 7.2	80 mM
sodium pyruvate	1.6 mM
NaCl	200 mM
NADH	0.20 mM
Stabilizers and preservatives	

STORAGE AND STABILITY

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The monoreagent is stable for 4 weeks at 2-8°C and for 1 week at room temperature ($\leq 25^\circ\text{C}$), when protected from the sunlight.

Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank ≤ 1.0

ADDITIONAL EQUIPMENT

General laboratory equipment.
Spectrophotometer; automated analyzer or photometer with thermostated cuvette, 1cm light path

SAMPLE

Serum, heparinized plasma. Samples free from hemolysis should be used. The enzyme in serum is stable for 2 days at 2-8°C. Frozen samples are rapidly inactivated.

CAUTION

The reagent contains sodium azide at 0.09%. Handle with care.
The safety statements are on the label. It is advisable to look at the SDS before using the reagent.
The disposal of the residues has to be made according to local regulations.

QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

PROCEDURE

The method described is the one proposed by the French Society of Clinical Biology. Bring the working reagent and the instrument to the working temperature. (30° / 37°C)

Monoreagent technique	30° / 37°C
Working Reagent	1.0 mL
Sample	0.02 mL

Bireagent technique	30° / 37°C
Buffer solution (A)	1.0 mL
Sample	0.02 mL
Mix, incubate for approx. 1 min	
NADH (B)	0.25 mL

Mix, read the absorbance after 1 min and start the stopwatch.
Read again the absorbance after 1, 2 and 3 min.

Determine the average value of $\Delta\text{Abs} / \text{min}$

Reading

Wavelength: 334 nm; 340 nm; 365 nm
Blank: Water
Cuvette: Thermostated; 1cm light-path

CALCULATIONS

The formula indicated is used to obtain the factor to calculate the U/L

$$\Delta\text{Abs}/\text{min} \times \frac{Vt \times 10^6}{\epsilon \times l \times Vs} = \text{U/L}$$

Where:

Vt: Total volume of reaction mixture;

Vs: Sample volume

l: Cuvette light path

ϵ : Extinction coefficient of NADH:

365 nm: 3.40×10^3

340 nm: 6.31×10^3

334 nm: 6.17×10^3

Determine the average value of $\Delta\text{Abs} / \text{min}$.

$$\text{U/L} = \Delta\text{Abs}/\text{min} \times \text{Factor}$$

Factors

	Monoreagent	Bireagent
334 nm	8252	10275
340 nm	8095	10080
365 nm	15000	18675

REFERENCE VALUES

Adults (37°C) 200 - 400 U/L

The stated values are for guidance. Each particular laboratory should establish its own normal range, using its own instrumentation, blood collection methods and test procedures

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using a manual method at 37°C and 340 nm.

Sensitivity, as detection limit: 10 U/L

Linearity: Up to 1800 U/L. For higher values, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.

Accuracy: 98.8 %

Repeatability, as Coefficient of Variation: 1.90%

Reproducibility, as Coefficient of Variation: 2.14%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

INTERFERENCES

It is recommended to use heparin as anticoagulant. Others, like citrate, oxalate, or fluoride may interfere with the test. Hemolyzed samples will give false results.

When assaying high activity samples, a very low initial absorbance can be found, which is mainly due to the rapid consumption of the NADH at the early stage of the reaction. In such a case, dilute the sample with saline (NaCl 0.9 %) and assay it once again.

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