COLESTEROL - LDL

MÉTODO CON SULFATO DE POLIVINILO

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Methods of Enzymatic Analysis, 3rd Edition, Volume VIII, Edited by H. U. Bergmeyer. (1985), 154 - 160.

Assmann, G., Jabs, H. U., Kohnert, U., Nolte, W., Schriewer, H., (1984), Clin. Chim. Acta., 140, 77 - 83.

Demacker, P. N., et al. (1984), Clin. Chem., 30, 1797 -1800.

Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

Zapico E. y Ordóñez J. (2002), Clin Invest Ar ter ioscl;14(5):272-6.

LDL - CHOLESTEROL

POLYVINYL SULPHATE METHOD

For "in vitro" determination of LDL - cholesterol in serum

PRINCIPLE

LDL-cholesterol can be determined as the difference between total cholesterol and the cholesterol content of the supernatant after precipitation of the LDL fraction by polyvinyl sulphate (PVS) in the presence of polyethylene-glycol monomethyl ether.

DIAGNOSTIC USE

LDL proteins are low density lipoproteins that transport cholesterol, triglycerides and lipids in the blood. A high level of cholesterol bound to the LDL fraction, LDL cholesterol, is associated with myocardial infarction and cardiovascular disease.

Other diseases that present high levels of LDL cholesterol are: nephrosis, diabetes, hypothyroidism or obesity.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

REAGENTS

1 x 10 mL Precipitant solution Ref. 99 06 10
Dropper for 100 tests

WORKING REAGENT PREPARATION

Reagent is ready to use.

WORKING REAGENT COMPOSITION

Concentrations in the reagent solution are:

Polyvinyl sulphate	0.7 g/L
EDTA Na ₂	5.0 mM
Polyethyleneglycol monomethyl ether	170 g/L
Stabilizers	

STORAGE AND STABILITY

When kept at 2 - 8°C, the reagent will remain stable until the expiration date stated on the label.

Signs of reagent deterioration:

Presence of particles. Non-homogeneous appearance.

ADDITIONAL EQUIPMENT

General laboratory equipment.
Spectrophotometer or photometer thermostable at 37°C.
Centrifuge.

Reagent for cholesterol determination and standard

It is advisable to use QCA reagents:

Kit 1 x 100 mL 99 52 82
Kit 3 x 100 mL 99 52 80
Kit 2 x 250 mL 99 50 12

Standard is included in the kits

CAUTION

The safety statements are on the label. It is recommended to read the MSDS before the reagent handling.
Waste products must be handled as per local regulations.

SAMPLE

Serum. The serum may be stored at 2 - 8°C for 1 day. Do not freeze.

QUALITY CONTROL

We recommend the inclusion of control sera Seriscann Normal (Ref. 99 41 48) and Seriscann Abnormal (Ref. 99 46 85) in each cholesterol determination to verify the results. For the precipitation reaction control, it is recommended to use liquid or lyophilized serum in which the native state of lipoproteins is retained. We recommend using a "pool" of normal sera of known value in LDL-cholesterol.
Each particular laboratory should establish its own control program.

AUTOANALYZERS

The precipitation reaction cannot be automated. Adaptations for cholesterol determination to different autoanalyzers are available on request.

INTERFERENCES

Sera with Triglyceride higher than 400 mg/dL could interfere. In such a case the sample shall be diluted 1/2 with saline (NaCl 0.9 %) prior to precipitation. Multiply the final result by 2.

Neither aged nor haemolysed sera will be assayed

PROCEDURE

1. Precipitation reaction

Precipitant solution 0.1mL (3 drops)
Sample 0.2 mL
Mix and let stand for 15 min at room temperature (20 - 25 °C).
Centrifuge at 2000 x g (1500 -2300 rpm)/ 15 min or
Centrifuge at 10000 x g (8000 -12000 rpm)/ 2 min
Determine cholesterol concentration in the supernatant.

2.- Determination of cholesterol

Samples: Total serum (SA) and supernatant (SP)

Technique	BL mL	SA mL	SP mL	ST mL
Serum	--	0.01	0.01	--
Supernatant	--	--	--	--
Standard	--	--	--	0.01
Working reagent	1.00	1.00	1.00	1.00

Mix well and let stand for 5 min at 37 °C or 10 min at room temperature.

Reading

Wavelength: 546 nm; 505 nm

Blank: the contents of BL

Colour stability: 1 hour(when protected from direct light)

CALCULATIONS

A. Cholesterol in supernatant

$$\frac{\text{Abs SA}}{\text{Abs ST}} \times 200 \times 1.5 = \text{mg supernatant cholesterol / dL}$$

B. Cholesterol in total serum

$$\frac{\text{Abs SP}}{\text{Abs ST}} \times 200 = \text{mg total cholesterol / dL}$$

Where:

Abs SP: Serum sample Absortion

Abs SA: Supernatant Absortion

Abs ST: Standard Absortion

C. LDL- cholesterol

$$\text{LDL- cholesterol (mg/dL)} = \text{total cholesterol in the sample (mg/dL)} - \text{supernatant cholesterol (mg/dL)}$$

SI Units: mg/dL x 0.0259 = mmol/L

Note: When the precipitant reagent is added to the sample, the latter is diluted by a factor of 1.5. Therefore, the cholesterol value should be multiplied by such a factor.

REFERENCE VALUES

Coronary heart disease risk:

< 100 mg/dL	Optimal
130 - 159 mg/dL	Borderline high
160 - 189 mg/dL	High
> 189 mg/dL	Very high

PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the precipitant reagent, the reagent to determine cholesterol and the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained by manual method. The cholesterol has been determinate with the reagent cholesterol liquid from QCA.

Sensitivity, as detection limit: 6 mg/dL

Accuracy: 96.4%

Repetitivity, as CV%: 2.87%

Reproducibility, as CV%: 2.96%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the comparison experiments are available on request.

REFERENCES

Methods of Enzymatic Analysis, 3rd Edition, Volume VIII, Edited by H. U. Bergmeyer. (1985), 154 - 160.

Assmann, G., Jabs, H. U., Kohnert, U., Nolte, W., Schriewer, H., (1984), Clin. Chim. Acta., 140, 77 - 83.

Demacker, P. N., et al. (1984), Clin. Chem., 30, 1797 -1800.

Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).

National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III), (2001) NIH Publication. Bethesda: National Heart, Lung, and Blood Institute.

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